

A rapid and simple method for screening fungi for extracellular protease enzymes

Thirunavukkarasu N¹, Suryanarayanan TS^{2*} , Rajamani $T^{1,2}$ and Paranetharan $MS^{1,2}$

¹ PG & Research Department of Botany, Ramakrishna Mission Vivekananda College, Chennai 600 004

² Vivekananda Institute of Tropical Mycology, Ramakrishna Mission Vidyapith, Chennai 600 004 *Corresponding authors t. sum 2002@values.com

*Corresponding author: t_sury2002@yahoo.com

Thirunavukkarasu N, Suryanarayanan TS, Rajamani T, Paranetharan MS 2017 – A rapid and simple method for screening fungi for extracellular protease. Mycosphere 8(1), 131–136, Doi 10.5943/mycosphere/8/1/12

Abstract

A simple and rapid method to detect extracellular protease produced by fungi is described. This method is more sensitive than the conventional agar plate method and could be used for screening a large number of samples in a short time. It also enables the visualization of different types of proteases. This method could be used with suitable modifications to detect other extracellular enzymes.

Keywords– Agar plate assay– enzyme assay – fungal enzymes – fungal protease

Introduction

Fungal proteases are used in the food, dairy, detergent, leather and pharmaceutical industries. They are also used for bioremediation and production of therapeutic peptides (Wu et al. 2006, Yike 2011). Fungal enzymes are more suited for industrial applications as fungi are easily cultured and are fast growing; many fungal enzymes are extracellular and are amenable to manipulations. Furthermore, fungal proteases are more diverse and exhibit wider substrate specificity (Kudryavtseva et al. 2008). Apart from engineering known enzymes for better efficiency by methods including directed evolution and rational design, screening microbes from relatively unexplored ecological niches for novel enzymes with desirable properties is a common method for improving the current industrial enzyme scenario (Kumar et al. 2011, Li et al. 2012). Fungal endophytes are non-disease causing endosymbionts of plants and comprise a group of different lineages (Guo et al. 2001, Hyde & Soytong 2008). Numerous studies have denoted the potential of endophytes to produce secondary metabolites exhibiting different bioactivities (Suryanarayanan et al. 2009, Kaushik et al. 2014, Mohana Kumara et al. 2014, Chen et al. 2016). Enzymes are produced in nature by endophytic fungi (Promputtha et al. 2010), yet they have not been explored intensively for their enzyme range (Suryanarayanan et al. 2012). Using endophytes, we describe here a simple and rapid qualitative assay for secreted fungal proteases.

Materials & Methods Fungal strains

Fungal endophytes residing in the inside the root, bark, or leaf tissues of some plants including forest trees (Western Ghats, Tamilnadu state, southern India) and mangrove forest (Pichavaram, Tamilnadu state, southern India) and maintained in the culture collection of VINSTROM were used in this study. Sixty one endophyte strains were screened to detect secreted protease enzymes by a qualitative dot blot method. This method was earlier used by us for detecting extracellular chitinase and chitosanase of endophytes (Govinda Rajulu et al. 2011); it was suitably modified for detecting proteases in the current study.

Detection of Protease enzyme

Preparation of samples

Each fungal isolate was grown in 100 ml Potato Dextrose medium (200 g Potato, 20 g Dextrose, 1000ml distilled water, pH 6.0) for five days as static culture at 26 °C. After incubation, the mycelium was filtered and the culture filtrate was dialyzed for 15 h against distilled water. The dialyzed culture filtrate was lyophilized and used as crude enzyme source for detecting proteases.

Gel preparation for spot assay and visualization of enzyme activity

A dot blot method described by Govinda Rajulu et al (2011) for detecting chitin modifying enzymes of fungi was modified suitably and used. A gel was prepared by mixing a solution containing 1 ml of 30 % acrylamide/bisacrylamide, 0.4 ml substrate (2 % autoclaved bacteriological gelatin – Qualigens 24345), 2.3 ml of 50mM TrisHCl buffer (pH 9.0), 0.003 ml of 100 % TEMED and 0.003 ml of 40 % ammonium persulphate and poured in a gel cassette. After polymerization, it was topped with the above mixture having Bis Tris buffer (pH 7.0) or sodium acetate buffer (pH 5.0). Thus, a composite gel made of gel strips of pH 5.0, pH 7.0 and pH 9.0 was obtained.10 mg of the lyophilized culture filtrate of each fungus (crude enzyme source) was mixed with 1 ml of appropriate buffer (buffers of different pH used for gel preparation) and centrifuged at 14000 rpm for 5 minutes (20 °C); 5 µl of the supernatant of this solution was spotted on the appropriate gel fields of the composite gel. The gel was incubated for 8-10 hours, stained with 0.025 % Coomassie Brilliant Blue (R 250) for 3 hours, washed with distilled water and then observed for the presence of clear zones on the deep blue gels which indicated enzyme activity. As controls, alkaline protease (Sigma P8038) and acidic protease (Sigma P0107) were dissolved in appropriate buffers and 5 µl of this solution (fresh or autoclaved) was spotted along with the samples on the gels.

Results

Protease action was clearly visible as a colourless spot on the dark blue background of the gel (Fig. 1). We observed that out of the 61 endophyte isolates, 25, 24 and 20 elaborated acidic, neutral and alkaline protease(s) respectively (Table 1). Interestingly, 17 isolates were positive for all the three types of proteases.

Discussion

Agar plate assays are most commonly used for detecting secreted proteases of fungi. This involves amending agar culture medium with a protein substrate such as casein or milk and growing the test fungus on this medium. Production of protease is visualized by the appearance of a clear zone around the fungal colony indicative of the hydrolysis of protein. Flooding the agar plate with a solution of Tri Chloro Acetic acid (Larsen et al. 1998) or Tannic acid (Saran et al. 2007) enhances the visibility of the cleared zone. In the dot blot assay, we used Coomassie Brilliant Blue stain which is routinely used for staining proteins separated by SDS-PAGE. This stain forms a complex with the protein and is sensitive enough to visualize proteins as low as 30 ng (Simpson 2010). It stains the gelatin substrate present in the gel deep blue and the clearing of this colour in the gel indicates protease action. The dot blot method is superior to the conventional agar plate assay in several ways. It is more sensitive and could be used for screening a large number of samples in a short time. By altering the pH of the supporting gel, it effectively detects different

AlternariaElacodendron glaucum / LAlternaria sp.Cordia obliqua / LAlternaria sp.Acacia chundra / LAureobasidium pullulansSpinifex littoreus / LAureobasidium pullulansSpinifex littoreus / LAureobasidium pullulansSpinifex littoreus / LChaetomium sp.Bruguiera cylindrica / LChaetomium sp.Bruguiera cylindrica / LChaetomium sp.Windentified / S-++Colletotrichum acutatumNeolitsea zeylanica / L+++Colletotrichum sp.Murraya paniculata /LColletotrichum sp.Kydia calycina / LColletotrichum sp.Rhizophora apiculata / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Conaga peltuta / LColletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LColletotri	Endophyte	Host	pH		
Alternaria sp.Cordia obliqua / LAlternaria sp.Acacia chundra / LAureobasidum pullulansSpinife's littoreus / LAureobasidum pullulansSpinife's littoreus / LBeltrania rhombicaSympiocos obusa / L+Chaetonium fbripilumGmeltna asiatica / BChaetonium sp.Bruguiera cylindrica / L-++Chaetonium sp.Nvicennia marina / R+++Collectorichum sp.Nyzgium cumini / LCollectorichum sp.Murraya paniculata /L+Collectorichum sp.Muraya paniculata /L+++Collectorichum sp.Smige anacardium / L+++Collectorichum sp.Rhizophora apiculata /L+++Collectorichum sp.Rhizophora apiculata /LCorynespora cassicolaButea monosperma /L+++Curvularia lunataIlex wightiana /BCurvularia sp.Acacia leucophloea /L+++Cylindrocladium sp.Vinca rosea /LCylindrocladium sp.Vinca rosea /LCylindrocladium sp.Vinca rosea /LCylindrocladium sp.Vinca rosea /LCylindrocladium sp.Vinca rosea /L- <th></th> <th></th> <th>5.0</th> <th></th> <th>9.0</th>			5.0		9.0
Alternaria sp.Cordia obliqua / LAlternaria sp.Acacia chundra / LAureobasidium pullulansSpinifes' litoreus / LAureobasidium pullulansSymplocos obusa / L+Beltrania rhombicaSymplocos obusa / L+Chaetomium sp.Bruguiera cylindrica / LChaetomium sp.Bruguiera cylindrica / LChaetomium sp.Unidentified / S-++Collectorichum sp.Syzygium cumin / LCollectorichum sp.Murraya paniculata ICollectorichum sp.Kydia calycina / L+++Collectorichum sp.Semecarpus anacardium / L+++Collectorichum sp.Rhizophora apiculata / LCollectorichum sp.Rhizophora apiculata / LCorynespora cassicolaButea monosperma / L+++Curvularia sp.Acacia leucophloea / LCyrularia sp.Nicaranga peltata / LCyrularia sp.Vinca rosea / LCyrularia sp.Cinnamonum malabatrum / L+++Cylindrocladium sp.Vinca rosea / LCylindrocladium sp.Vinca rosea / LCylindrocladium sp.Vinca rosea / L- </td <td>Alternaria alternata</td> <td>Elaeodendron glaucum / L</td> <td>-</td> <td>-</td> <td>-</td>	Alternaria alternata	Elaeodendron glaucum / L	-	-	-
Alternaria sp.Acacia chundra / LAureobasidium pullulansSpinifex littoreus / LBeltrania rhombicaSymplocos obtusa / L+Beltrania rhombicaSymplocos obtusa / L+Chaetomium fibripiliumGmelina asiatica / BChaetomium sp.Bruguiera cylindrica / LChaetomium sp.Wticennia marina / R+++Collectorichum acutatumNeolitsea czylanica / LColletorichum sp.Murraya paniculata /LColletorichum sp.Murraya paniculata /LColletorichum sp.Murraya paniculata /L+++Colletorichum sp.Rhizophora apiculata /L+++Colletorichum sp.Rhizophora apiculata /L+++Colletorichum sp.Rhizophora apiculata /L+++Colletorichum sp.Rhizophora apiculata /L+++Colletorichum sp.Syzygium calophyllifolium /LCurvularia lunataIlex wightiana / BCurvularia sp.Macaranga peltata /L++++Cyridorcladium sp.Cinnamomum malabatrum /L++++Fusarium sp.Pongamia pinnata /BFusarium sp.Pongamia pinnata /B	<i>Alternaria</i> sp.	0	-	-	-
Aureobasidium pullulansSpinifex littoreus / LA. pullulansEuonymus angulatus / LA. pullulansEuonymus angulatus / LChaetomium fibripiliumGmelina asiatica / BChaetomium sp.Bruguiera cylindrica / LChaetomium sp.Bruguiera cylindrica / LChaetomium sp.Avicennia marina / R++Collectorichum sp.Muraya paniculata / LCollectorichum sp.Muraya paniculata / LCollectorichum sp.Kydia calycina / LCollectorichum sp.Kydia calycina / LCollectorichum sp.Rhizophora apiculata / LCollectorichum sp.Rhizophora apiculata / LCollectorichum sp.Rhizophora apiculata / LCorynespora cassicolaTerminalia crenulata / LCryptosporiopsis sp.Syzyjum calophyllfolium / LCurvularia lunataIlex wightinan / BCurvularia sp.Macaica leucophloea / L++Fusarium sp.Pongamia pinnata / BOrechslera sp.Vinca rosca / LOrechslera sp.Vinca rosca / LCylindrocladium sp.Glochidin sp. / B++Histophora apiculata / LChaetorichum sp.Pongamia pinnata / BChaetoric	1	1	-	-	-
A. pullulansÉuonymus angulatus / LBeltrania rhombicaSymplocos obtusa / L+Chaetomium sp.Bruguiera cylindrica / LChaetomium sp.Unidentified / S-++Chaetomium sp.Unidentified / S-++Collectorichum acutatumNeolitsea zeylanica / L+++Collectorichum sp.Syzygium cumini /LCollectorichum sp.Murraya paniculata /LCollectorichum sp.Murraya paniculata /LCollectorichum sp.Rhizophora apiculata /LCollectorichum sp.Rhizophora apiculata /LCollectorichum sp.Rhizophora apiculata /LCorynespora cassicolaButea monosperma /L+++Cryptosporiopsis sp.Syzygium calophyllifolium /LCurvularia sp.Acacia leucophloea /L+++Cylindrocladium sp.Cinnamonum malabatrum /L+++Fusarium sp.Pongamia pinnata /BGolmerella cingulata cingulata sonnifera /LCurvularia sp.Acacia leucophloea /L++++Corgonspors oryzaeUnidentified /LSemium sp.Philophora apiculata /L <tr <tr="">Cyl</tr>	1	Spinifex littoreus / L	-	-	-
Beltrania rhombicaSymplocos obtusa / L+Chaetomium fibripiliumGmelina asiatica / BChaetomium sp.Bruguiera cylindrica / LChaetomium sp.Avicennia marina / R+++Collectorichum acutatumNeolitsea zeylanica / L+++Collectorichum sp.Syzygium cumini / LCollectorichum sp.Murraya paniculata / LCollectorichum sp.Kydia calycina / LCollectorichum sp.Rhizophora apiculata / L+++Collectorichum sp.Rhizophora apiculata / LCollectorichum sp.Rhizophora apiculata / LCorynespora cassicolaTerminalia crenulata / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / LCurvularia sp.Cinnamonum malabatrum / L+++-Fusarium sp.Pongamia pinnata / BFusarium sp.Rhizophora apiculata / LCurvularia sp.Cinnamonum malabatrum / L+++-Curvularia sp.Cinnamonum malabatrum / L+++-Fusarium sp.Pongamia pinnata / B	-		-	-	-
	-		+	-	-
Chaetomium sp.Bruguiera cylindrica / LChaetomium sp.Unidentified / S-++Chaetomium sp.Avicennia marina / R+++Colletotrichum sp.Nycennia marina / R+++Colletotrichum sp.Murraya paniculata / LColletotrichum sp.Murraya paniculata / LColletotrichum sp.Unidentified / L+++Colletotrichum sp.Semecarpus anacardium / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Corprespora cassicolaTerminalia crenulata / LC. cassicolaButea monsperma / L++++Curvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L++++Cylindroctadium sp.Cinnamonum malabatrum / L++-Drechslera sp.Vinca rosea / LClimarium sp.Pongamia pinnata / BOdilisporium sp.Glochidin sp. / B++++NorgaeDaphiphyllum neilgherrense / LNodulisporium sp.Glochidin sp. / B++++Phoma sp.Avicennia marina / R++++Phomops			-	-	-
Chaetomium sp.Unidentified / S-++Chaetomium sp.Avicennia marina / R+++Colletotrichum acutatumNeolitsea zeylanica / L+++Colletotrichum sp.Syzgium cumini /LColletotrichum sp.Murraya paniculata /LColletotrichum sp.Kydia calycina / LColletotrichum sp.Unidentified / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Colletotrichum sp.Rhizophora apiculata / LCornespora cassicolaTerminalia crenulata / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Cylindrocladium sp.Cinnamonum malabatrum / L+++Fusarium sp.Pongamia pinnata / BFusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nodulisporium sp.Philzophora apiculata / R+++NograeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidon sp. / B++++Phomogsis	• •		-	-	-
Chaetomium sp.Avicennia marina / R+++Colletotrichum acutatumNeolitsea zeylanica / L+++Colletotrichum sp.Syzygium cumini / LColletotrichum sp.Murraya paniculata / L+Colletotrichum sp.Unidentified / L+++Colletotrichum sp.Unidentified / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / LCorynespora cassicolaTerminalia crenulata / LCurvularia lantaIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Acacia leucophloea / L+++Cylindrocladium sp.Cinnamomum malabatrum / L++-Drechslera sp.Vinca rosea / LNigrospora oryzaeUnidentified / L+++Nodulisporium sp.Pongamia pinnata / BNodulisporium sp.Glochidin sp. / B+++Prestalotiopsis sp.Kitaphora apiculata / LNodulisporium sp.Glochidin sp. / B+++Phoma sp.Phogamia pinnata / BPhoma sp.Vincania sannifera / L <td< td=""><td>÷</td><td>e ,</td><td>-</td><td>+</td><td>-</td></td<>	÷	e ,	-	+	-
Colletotrichum acutatumNeolitsea zeylanica / L+++Colletotrichum sp.Syzygium cumini /LColletotrichum sp.Murraya paniculata /LColletotrichum sp.Unidentified / L+++Colletotrichum sp.Semecarpus anacardium / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Corprespora cassicolaTerminalia crenulata / LCryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L++++Cylindrocladium sp.Mina pontanta / BDrechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata / R++++NoryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidion sp. / B+++Nodulisporium sp.Glochidion sp. / B+++Phoma sp.Avicennia marina / R+++Phoma sp.Avicennia marina / R+++ </td <td>-</td> <td></td> <td>+</td> <td></td> <td>+</td>	-		+		+
Colletotrichum sp.Syzygium cumini /LColletotrichum sp.Murraya paniculata /L+Colletotrichum sp.Kydia calycina /L+++Colletotrichum sp.Semecarpus anacardium /L+++Colletotrichum sp.Rhizophora apiculata /L+++Colletotrichum sp.Rhizophora apiculata /LColletotrichum sp.Rhizophora apiculata /L+++Corprespora cassicolaButea monosperma /L+++Cryptosporiopsis sp.Syzygium calophyllifolium /LCurvularia lunataIlex wightiana /BCurvularia sp.Acacia leucophloea /L+++Cylidrocladium sp.Cinnamonum malabatrum /L++-Drehslera sp.Vinca rosea /LFusarium sp.Pongamia pinnata /BGlomerella cingulataCanthium parviflorum /L+++NoryzeDaphinphyllum neilgherrense /LNodulisporium sp.Glochidion sp. /B++++Pendulisporium sp.Phyllonta aficalata /L+++Phomopsis sp.Acacia suma/LNodulisporium sp.Glochidin sp. /B++++Phomopsis sp.Acacia suma/LPhoma sp.Unid	-		+	+	+
Colletotrichum sp.Murraya paniculata /L+Colletotrichum sp.Kydia calycina / LColletotrichum sp.Unidentified / L++++Colletotrichum sp.Rhizophora apiculata / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Corynespora cassicolaTerminalia crenulata / L+++Corynespora cassicolaTerminalia crenulata / LC cassicolaButea monosperma / L+++Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia sp.Acacia leucophloea / L+++Cylindrocladium sp.Cinnamomum malabatrum / L++-Cylindrocladium sp.Cinnamomum malabatrum / LFusarium sp.Pongania pinnata / BGlomerella cingulataCanthium parviflorum / L+++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidion sp. / B+++Pestalotiopsis sp.Avicennia marina / R+++Phomoa sp.Avicennia marina / R+++Phomos sp.Avicennia marina / R+++Phomopsis sp.Avicennia marina / R+++ <trr></trr>		-	_	_	_
Colletotrichum sp.Kydia calycina / LColletotrichum sp.Unidentified / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Corynespora cassicolaTerminalia crenulata / LC. cassicolaButea monosperma / L++++Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L++++Cylindrocladium sp.Cinnamomum malabatrum / L+++-Drechslera sp.Vinca rosea / LFusarium sp.Pongamia pinnata / BNoidlisporium sp.Pongamia pinnata / BNodulisporium sp.Qlachidium sp.Glochidions p. / B+++Nodulisporium sp.Glochidions p. / B++++Nodulisporium sp.Glochidions p. / BPhoma sp.Avicennia marina / R+++++Phoma sp.Avicennia marina / R++++Phomopsis sp.Avicennia marina / R++++Phomopsis sp.Avicennia marina / R++++ <tr <tr="">Phomopsis</tr>	-		_	_	+
Colletotrichum sp.Unidentified / L+++Colletotrichum sp.Semecarpus anacardium / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Cornespora cassicolaTerminalia crenulata / LC. cassicolaButea monosperma / L+++Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophoea / L+++Curvularia sp.Macaranga peltata / L-+-Cylindrocladium sp.Cinnamomum malabatrum / L+++Fusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongania pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nodulisporium sp.Glochidion sp. / B+++Nodulisporium sp.Glochidion sp. / B+++Pencillium sp.Phyllanthus emblica / BNodulisporium sp.Rhizophora apiculata / LNodulisporium sp.Glochidion sp. / B+++Pencillium sp.Phyllanthus emblica / BPhoma sp.Avicennia marina / R+++Phomopsis sp.Avicennia marin	±		_	_	_
Colletotrichum sp.Semecarpus anacardium / L+++Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Corynespora cassicolaTerminalia crenulata / LC. cassicolaButea monosperma / L+++Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / L-+-Curvularia sp.Cinnamonum malabatrum / L++-Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata / BGlomerella cingulataCanthium parviflorum / L+++NograpeDaphniphyllum neilgherense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Pongania pinnata / R++++Nodulisporium sp.Glochidion sp. / BNodulisporium sp.Unidentified / LPostalotiopsis sp.Avicennia marina / R++++Phoma sp.Av	-	• •	+	+	+
Colletotrichum sp.Rhizophora apiculata / LColletotrichum sp.Rhizophora apiculata / L+++Corynespora cassicolaTerminalia crenulata / L++Crynespora cassicolaButea monosperma / L++Crynosporiopis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L++Curvularia sp.Macaranga peltata / L-+Cylindrocladium sp.Cinnamomum malabatrum / L+-Drechslera sp.Vinca rosea / LFusarium sp.Pongamia pinnata / BFusarium sp.Pongamia pinnata / BNoryzaeUnidentified / L++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Withania somifera / LPestalotiopsis sp.Litsea stocksii / L++Phoma sp.Unidentified / S++Phoma sp.Unidentified / S+-Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R++Phomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / L<	1				
Colletotrichum sp.Rhizophora apiculata / L+++Corynespora cassicolaTerminalia crenulata / LC. cassicolaButea monosperma / L+++Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / L-+-Orechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidion sp. / B+++Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+++Phoma sp.Unidentified / S+Phoma sp.Unidentified / S+Phomopsis sp.Acacia suma / LPhomopsis sp.Avicennia marina / R++++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhomopsis sp.Avicennia marina / R++ <t< td=""><td>±</td><td>*</td><td>_</td><td>-</td><td>_</td></t<>	±	*	_	-	_
Corynespora cassicolaTerminalia crenulata / LC. cassicolaButea monosperma / L+++C. cassicolaButea monosperma / L+++Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / L-+-Curvularia sp.Cinnamomum malabatrum / L+Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++Nodulisporium sp.Glochidion sp. / B+++Pestalotiopsis sp.Rhizophora apiculata / LPestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Withania somifera / LPhoma sp.Rhizophora apiculata / L+++Phomopsis sp.Accia suma / L+++Phomopsis sp.Accia suma / LPhoma sp.Avicennia marina / R+++Phomopsis sp.Accia suma / LPhomops	-	· ·	+	+	+
C. cassicolaButea monosperma / L+++Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / L-++Curvularia sp.Macaranga peltata / L-+-Cylindrocladium sp.Cinnamonum malabatrum / L+++Cylindrocladium sp.Vinca rosea / LFusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++N oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidion sp. / B+++Pestalotiopsis sp.Rhizophora apiculata / L++Phoma sp.Phyllanthus emblica / BPhoma sp.Unidentified / S+++Phoma sp.Nicennia marina / R++Phomopsis sp.Avicennia marina / R++Phomopsis sp.Avicennia marina / R++Phomopsis sp.Acacia suma / LPhomopsis sp.Avicennia marina / R++Phomopsis sp.Avicennia marina / R++Phomopsis sp.Canthium parviflorum / LPithonyces gramincolaLitsea floribunda / LPithonyces g	-	· ·	-	-	-
Cryptosporiopsis sp.Syzygium calophyllifolium / LCurvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / L-+-Cylindrocladium sp.Cinnamonum malabatrum / L+Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++Nodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Nyllanthus emblica / BPhoma sp.Nicennia marina / R+++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LPithomyces graminicolaLitsea floribunda / LPithomyces graminicolaLitsea	· · ·		- _	-	- _
Curvularia lunataIlex wightiana / BCurvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / L-++Curvularia sp.Cinnamomum malabatrum / L+Cylindrocladium sp.Cinnamomum malabatrum / L+Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata / R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidion sp. / B+++Pestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Phyllanthus emblica / BPhoma sp.Unidentified / S+++Phoma sp.Unidentified / S+++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LP. capitalensisBridelia retusa / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LP. capitalensisOugeinia oojeinensis / LP. capitalensis <td< td=""><td></td><td>-</td><td>Т</td><td>Т</td><td>Т</td></td<>		-	Т	Т	Т
Curvularia sp.Acacia leucophloea / L+++Curvularia sp.Macaranga peltata / L-+-Cylindrocladium sp.Cinnamomum malabatrum / L+Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BNodulisporium sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Navicennia marina / R+++Phomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LP. capitalensisBridelia retusa / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++TanialensisIltsea floribunda / L			-	-	-
Curvularia sp.Macaranga peltata / L-+-Cylindrocladium sp.Cinnamomum malabatrum / L+Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B+++Pericillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Unidentified / S+++Phoma sp.Avicennia marina / R+++Phomopsis sp.Accaia suma / LPhomopsis sp.Acaia suma / LPhomopsis sp.Acaia floribunda / LP. capitalensisBridelia retusa / LP. capitalensisDugenia oojeiensis / LP. capitalensisDugenia oojeiensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++ <tr tr=""><</tr>		e e	-	-	-
Cylindrocladium sp.Cinnamomum malabatrum / L+Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++NoryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Glochidion sp. / B+++Pericillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Nitera stocksii / L+++Phoma sp.Unidentified / S+Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LP. capitalensisBridelia retusa / LP. capitalensisDugenia oojeinensis / LP. capitalensisDugenia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++Hizophora apiculata / R++++Pitelie formRhizophora apiculata / R+++Hizophora apiculat	-	•	+		+
Drechslera sp.Vinca rosea / LFusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++N. oryzaeDaphniphyllum neilgherense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Unidentified / S+++Phoma sp.Unidentified / S+Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhomopsis sp.Canthium parviflorum / LPhomopsis sp.Acacia suma / LPhyllosticta capitalensisBridelia retusa / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia mar	-		-	Ŧ	-
Fusarium sp.Rhizophora apiculata/ R+++Fusarium sp.Pongamia pinnata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Unidentified / S+Phoma sp.Unidentified / S+Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++Piterile formRhizophora apiculata / R+++Hex wightiana / LPiterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / LIter formRhizophora apiculata / R+++Iter f	· ·		+	-	-
Fusarium sp.Pongamia pinata / BGlomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L+++N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Phoma sp.Unidentified / S+Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhomopsis sp.Bridelia retusa / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R++++Sterile formRhizophora apiculata / R++++Taeniolella sp.Ilex wightiana / LSterile formRhizophora apiculata / R	-		-	-	-
Glomerella cingulataCanthium parviflorum / L+++Nigrospora oryzaeUnidentified / L++-N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++Hitting application opieinensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++Sterile formAvicennia marina / R+++Taeniolella sp.Ilex wightiana / LSterile formAvicennia marina / R+++Taeniolella sp.Ilex wightiana / L<	I I		+	+	+
Nigrospora oryzaeUnidentified / L++-N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisCanthium parviflorum / LSterile formAvicennia marina / R++++Hyllosticta sp.Litsea floribunda / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++++Sterile formRhizophora apiculata / R++++Taeniolella sp.Ilex wightiana / LSterile formRhizophora apiculata / R+++++Paniolella sp.Ilex wightiana / L- <td>-</td> <td>· ·</td> <td>-</td> <td>-</td> <td>-</td>	-	· ·	-	-	-
N. oryzaeDaphniphyllum neilgherrense / LNodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L++Pestalotiopsis sp.Litsea stocksii / L++Phoma sp.Unidentified / S+-Phoma sp.Vnidentified / S+-Phoma sp.Avicennia marina / R++Phomopsis sp.Acacia suma / LPhomopsis sp.Acacia suma / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R++H++-Pitelia sp.Canthium parviflorum / LSterile formRhizophora apiculata / R++HSterile form <t< td=""><td>6</td><td>· ·</td><td></td><td></td><td>+</td></t<>	6	· ·			+
Nodulisporium sp.Withania somnifera / LNodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R+++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisCanthium parviflorum / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Hex wightiana / L	· ·		+	+	-
Nodulisporium sp.Glochidion sp. / B+++Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R+++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Hizophora apiculata / R+++Sterile formRhizophora apiculata / R++Taeniolella sp.Ilex wightiana / L	-		-	-	-
Penicillium sp.Phyllanthus emblica / BPestalotiopsis sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R+++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R++Taeniolella sp.Ilex wightiana / L		0	-	-	-
Pestalotiopsis sp.Rhizophora apiculata / L+++Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R+++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L		-	+	+	+
Pestalotiopsis sp.Litsea stocksii / L+++Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R++-Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LP. capitalensisCanthium parviflorum / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Taeniolella sp.Ilex wightiana / L	1	•	-	-	-
Phoma sp.Unidentified / S+Phoma sp.Avicennia marina / R+++Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LP. capitalensisOugeinia oojeinensis / LP. capitalensisCanthium parviflorum / LSterile formAvicennia marina / R+++Taeniolella sp.Ilex wightiana / L	· ·				
Phoma sp.Avicennia marina / R++-Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LP. capitalensisOugeinia oojeinensis / LP. capitalensisCanthium parviflorum / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Taeniolella sp.Ilex wightiana / L	· ·		+	+	+
Phomopsis sp.Avicennia marina / R+++Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Taeniolella sp.Ilex wightiana / L	1			-	-
Phomopsis sp.Acacia suma / LPhyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L	-		+	+	-
Phyllosticta capitalensisTurpinia nepalensis / LP. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Taeniolella sp.Ilex wightiana / L			+	+	+
P. capitalensisBridelia retusa / LP. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L			-	-	-
P. capitalensisOugeinia oojeinensis / LPithomyces graminicolaLitsea floribunda / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L			-	-	-
Pithomyces graminicolaLitsea floribunda / LSordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L	-		-	-	-
Sordaria sp.Canthium parviflorum / LSterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L	-		-	-	-
Sterile formAvicennia marina / R+++Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L			-	-	-
Sterile formRhizophora apiculata / R+++Taeniolella sp.Ilex wightiana / L	-	· ·	-	-	-
Taeniolella sp.Ilex wightiana / L			+	+	+
• •			+	+	+
Talaromyces sp.Avicennia marina / R+	-		-	-	-
	Talaromyces sp.	Avicennia marina / R	-	-	+

 Table 1 Detection of protease activity of endophytes by dot blot assay.

Endophyte	Host	pH		
		5.0	7.0	9.0
Torulomyces sp.	Randia dumetorum / L	-	-	-
Trichoderma sp.	Azadirachta indica / B	-	-	-
Trichoderma sp.	Pedilanthus tithymaloides / L	-	-	-
Trichoderma sp.	Unidentfied / L	-	-	-
Xylariaceous form	Cassia fistula / L	+	-	+
Xylariaceous form	Anogeissus latifolia / L	+	+	+
Xylariaceous form	Stereospermum tetragonum / L	-	-	-
Xylariaceous form	Syzygium calophyllifolium / L	-	-	-
Xylariaceous form	Ziziphus jujuba / L	+	+	+
Xylariaceous form	Butea monosperma/ L	+	+	-
Xylariaceous form	Dalbergia latifolia / L	-	+	-
Xylariaceous form	Semecarpus anacardium / B	+	+	-
Positive strains		25	24	20

L = isolated from leaf, B = isolated from bark, S = isolated from stem, R = isolated from root

Table 1 continued

proteases. Thus, large numbers of samples could be screened rapidly and potential producers couldbe identified for further detailed studies and use. Furthermore, the dot blot method could be modified suitably for detecting many other secreted enzymes such as amylase, and l-asparaginase. Recent studies show that fungal endophytes are a source of therapeutic and industrial enzymes such as chitin modifying enzymes (Govinda Rajulu et al. 2011), l-asparaginase (Nagarajan et al. 2014) and xylan-degrading enzymes (Thirunavukkarasu et al. 2015). Our present study shows that they could be a novel source of different proteases.




Figure 1 – Detection of protease activity of endophytes by dot blot assay. **a** Acidic protease sample. **b** Autoclaved acidic protease sample. **c**, **h**, **m** *Pestalotiopsis* sp. (*Litsea stocksii* / L). **d**, **i**, **n** *Pithomyces graminicola* (*Litsea floribunda* / L). **e**, **j**, **o** *Curvularia* sp. (*Acacia leucophloea* / L). **f**, **g** Not spotted. **k** Alkaline protease sample. **l** Autoclaved alkaline protease sample.

Acknowledgements

TSS thanks the Secretary, Ramakrishna Mission Vidyapith, Chennai, for providing research facilities. NT thanks the Department of Science and Technology, Government of India for funds through the Science and Engineering Research Board project (SB/EMEQ-005/2013).

References

- Chen C, Hu SY, Luo DQ, Zhu SY, Zhou CQ. 2016 Potential antitumor agent from the endophytic fungus *Pestalotiopsis photiniae* induces apoptosis via the mitochondrial pathway in HeLa cells. Oncology Reports 30, 1773–1781.
- Govinda Rajulu MB, Thirunavukkarasu N, Suryanarayanan TS, Ravishankar JP, El Gueddari NE, Moerschbacher BM. 2011 Chitinolytic enzymes from endophytic fungi. Fungal Diversity 47, 43–53.
- Guo LD, Hyde KD, Liew ECY 2001 Detection and taxonomic placement of endophytic fungi within frond tissues of *Livistona chinensis* based on rDNA sequences. Molecular Phylogenetics & Evolution 20, 1–13.
- Hyde KD, Soytong K 2008 The fungal endophyte dilemma. Fungal Diversity 33: 163-173.
- Kaushik NK, MuraliTS, Sahal D, Suryanarayanan TS. 2014 A search for antiplasmodial metabolites among fungal endophytes of terrestrial and marine plants of southern India. Acta Parasitologica 59, 745–757.
- Kudryavtseva OA, Dunaevsky YE, Kamzolkina OV, Belozersky MA. 2008– Fungal proteolytic enzymes: Features of the extracellular proteases of xylotrophic basidiomycetes. Microbiology 77, 643–653.
- Kumar L, Awasthi G, Singh B. 2011– Extremophiles: a novel source of industrially important enzymes. Biotechnology 10, 121–135.
- Larsen MD, Kristiansen KR, Hansan TK. 1998 Characterization of the proteolytic activity of starter cultures of *Penicillium roqueforti* for production of blue veined cheeses. International Journal of Food Microbiology 43, 215–221.
- Li S, Yang X, Yang S, Zhu M, Wang X. 2012 Technology prospecting on enzymes: application, marketing and engineering. Computational and Structural Biotechnology Journal 2, 1–11.
- Mohana Kumara P, Soujanya KN, Ravikanth G, Vasudeva R, Ganeshaiah KN, Uma Shaanker R. 2014 – Production of the chromane alkaloid, rohitukine and its attenuation in endophytic fungi isolated from *Dysoxylum binectariferum* Hook.f and *Amoora rohituka* (Roxb).Wight & Arn. Phytomedicine 21, 54 –546.
- Nagarajan A, Thirunavukkarasu N, Suryanarayanan TS, Gummadi SN. 2014 Screening and isolation of novel glutaminase free l-asparaginase from fungal endophytes. Research Journal of Microbiology 9, 163–176.
- Promputtha I, Hyde KD, McKenzie EHC, Peberdy JF; Lumyong, S. 2010 Can leaf degrading enzymes provide evidence that endophytic fungi becoming saprobes? Fungal Diversity 41, 89–99.
- Saran S, Isar J, Saxena RK. 2007 A modified method for the detection of microbial proteases on agar plates using Tannic acid. Journal of Biochemical and Biophysical Methods 70, 697–699.
- Simpson RJ. 2010 Rapid Coomassie blue staining of protein gels. Cold Spring Harbor protocols doi: 10.1101/pdb.prot5413.
- Suryanarayanan TS, Thirunavukkarasu N, Govindarajulu MB, Sasse F, Jansen R, Murali TS. 2009 – Fungal endophytes and bioprospecting. Fungal Biology Reviews 23, 9–19.
- Suryanarayanan TS, Thirunavukkarasu N, Govinda Rajulu MB, Gopalan V. 2012 Fungal endophytes: An untapped source of biocatalysts. Fungal Diversity 54, 19–30.
- Thirunavukkarasu N, Jahnes B, Broadstock A, Govinda Rajulu MB, Murali TS, Gopalan V, Suryanarayanan TS. 2015 Screening marine-derived endophytic fungi for xylan-degrading enzymes. Current Science109:112–120.

- Wu TY, Mohammada AW, Jahim JM. 2006 Investigations on protease production by a wild-type *Aspergillus terreus* strain using diluted retentate of pre-filtered palm oil mill effluent (POME) as substrate. Enzyme and Microbial Technology 39, 1223–1229.
- Yike I. 2011 Fungal proteases and their pathophysiological effects. Mycopathologia 171, 299–323.